

Guided Surgery

Introducing the adinguide™: a guided surgery tool designed to enhance your clinic.

Our kits are tailored to support effortless and precise implant placement. With our proprietary ActiveFlow™ Irrigation technology, where coolant is efficiently delivered to the bone, preventing overheating.

The adinguide™ drills feature carved channels for direct coolant flow from your hand-piece to the osteotomy site. Our intuitive drill kit accommodates different implant sizes and lengths, while the tissue punch and pilot drill allow flapless access.

The keyless procedure offers easy handling and clear visibility of the surgical site.

Our self-centering drills utilize a diameter cylinder at the top to ensure accurate positioning within the guide's sleeve. Additionally, built-in stoppers guarantee precise and consistent drilling depths.



Guided Surgery Kit

Implants Drills | GSK001



For fully guided procedures:




Unified Conical
Implant Platform



GS Tools
Kit UniFit
GSK005



5.5mmD Sleeve
GS0055



GS5055 Tissue Punch 5.0mmD, S-5.5



GS0825 Drill S-5.5, 2.5/1.8mmDx8.0mmL



GS1325 Drill S-5.5, 2.5/1.8mmDx13.0mmL



GS0625 Pilot Drill S-5.5, 2.5/1.8mmDx6.0mmL




GS1025 Drill S-5.5, 2.5/1.8mmDx10.0mmL



GS1625 Drill S-5.5, 2.5/1.8mmDx16.0mmL



GS1125 Drill S-5.5, 2.5/1.8mmDx11.5mmL



GS0632 Drill S-5.5, 3.2/2.5mmDx6.25mmL




GS0638 Drill S-5.5, 3.85/3.2mmDx6.25mmL



GS0646 Drill S-5.5, 4.6/3.85mmDx6.25mmL



GS0832 Drill S-5.5, 3.2/2.5mmDx8.0mmL



GS0838 Drill S-5.5, 3.85/3.2mmDx8.0mmL



GS0846 Drill S-5.5, 4.6/3.85mmDx8.0mmL



GS1032 Drill S-5.5, 3.2/2.5mmDx10.0mmL



GS1038 Drill S-5.5, 3.85/3.2mmDx10.0mmL



GS1046 Drill S-5.5, 4.6/3.85mmDx10.0mmL



GS1132 Drill S-5.5, 3.2/2.5mmDx11.5mmL



GS1138 Drill S-5.5, 3.85/3.2mmDx11.5mmL



GS1146 Drill S-5.5, 4.6/3.85mmDx11.5mmL



GS1332 Drill S-5.5, 3.2/2.5mmDx13.0mmL



GS1338 Drill S-5.5, 3.85/3.2mmDx13.0mmL



GS1346 Drill S-5.5, 4.6/3.85mmDx13.0mmL



GS1632 Drill S-5.5, 3.2/2.5mmDx16.0mmL



GS1638 Drill S-5.5, 3.85/3.2mmDx16.0mmL



GS1646 Drill S-5.5, 4.6/3.85mmDx16.0mmL

Guided Surgery Tools Kit UniFit Implants | GSK005



For fully guided procedures:




Unified Conical
Implant Platform



GS Drills Kit
GSK001


5.5mmD Sleeve
GS0055

1

GS1915 

Anchoring Drill S-2.0, 2.0/1.5mmDx19mmL

2

GS1917 

Anchoring Drill S-2.0, 2.0/1.7mmDx19mmL

3

GS0028 

Handpiece Guide Fixation Screw Driver

4

GS5041 

Guided UniFit Handpiece Fixture Mount S-5.5

5


GS5043 

Guided UniFit Ratchet Fixture Mount S-5.5

6

GS5037 

Guided UniFit Ratchet Fixture Mount S-5.5, Screw Retained (qty.4)



7

GS0023 

Manual Socket L10.5

9

GS0021 

Handpiece Socket L10.5

11

RP0004 

Star Hand Driver - Long

8

GS0024 

Manual Socket L14.5

10

GS0022 

Handpiece Socket L14.5

12

GS5046 

UniFit Ratchet Retrieval Tool

13

GS0027 

Guide Fixation Screw Driver (Surgical Driver)





14

GS6111 





Surgical Torque Ratchet 35-100 Ncm for hex 4mm

Components

Accessories

| | | | Catalog |
|---|--------------------------------------|-----------|-----------|
|  | Guided Vertical Fixation Screw S-5.5 | x 3 pack | GS0007 |
|  | 2.0mmD Anchoring Screw | x 3 pack | GS1920 |
|  | 5.5mmD Sleeve | x 1 pack | GS0055 |
| | | x 3 pack | GS0055x3 |
| | | x 10 pack | GS0055x10 |
|  | 2.0mmD Anchoring Sleeve | x 3 pack | GS0020 |

18mm Drills

| | | | |
|---|---|---------|--------|
|  | Drill S-5.5, 2.5/1.8mm Dx18.0mmL | Catalog | GS1825 |
|  | Drill S-5.5, 3.2/2.5mmD x18.0mmL | Catalog | GS1832 |
|  | Drill S-5.5, 3.85/3.2mmD x18.0mmL | Catalog | GS1838 |
|  | Drill S-5.5, 4.6/3.85mmD x18.0mmL | Catalog | GS1846 |